

Unisol W



Technical Data Sheet

Intended use: Additive for the improved wax-solving for the use in boiling-out units

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical designation
1	distilled water
2	sodium benzoate
3	fatty alcohol
4	terpene

Properties:

Colour: blue
Form: liquid
Odour: characteristic
Density: 1.0 – 1.1g/cm³
Boiling point: 100 °C

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 82909
dentamid@dreve.de
www.dreve.com