

Fixtemp cement



Technical Data Sheet

Product description: Eugenol-free, temporary cement

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemische Bezeichnung
1	zinc oxide
2	natural resin
3	fat acids
4	additives
5	paraffine/ vaseline

Processing properties:

Processing time at 23°C: \geq 1:00 min
Curing time in the mouth: approx. 4:00 – 7:00 min
Removal of surplus material
after setting: approx. 5:00 – 6:00 min

Physical properties:

Compressive strength: \leq 35 MPa
Film- thickness: \leq 20 μ m
Acid- soluble as per content: \leq 1 mg/kg
Decomposition (after 24h, m/m): \leq 2.5 %

Conforms to ISO 3107: 2011 Dentistry:

Zinc oxide / eugenol cements and zinc oxide / non-eugenol cements, Type1, Class 1

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com/dentamid