

# Regofix opaque



## Technical Data Sheet

**Intended use:** Bite registration material on vinylpolysiloxane basis.

**Manufacturer:** Drewe Dentamid GmbH  
Max-Planck-Straße 31  
59423 Unna / Germany

### Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Polydimethylsiloxane Oil
3	Filling materials
4	Platinum catalyst
5	Polymethylhydrogensiloxane
6	Pigments

### Properties:

Colour: blue

Hardness Shore A:  $90 \pm 5$

Consistency: heavy - body - Type 1 (ISO 4823)

Elastic recovery:  $> 99.8 \%$  (ISO 4823)

Linear dimensional change:  $< 0.2 \%$  (ISO 4823)

  

Processing time:  $30 \pm 5$  sec at  $23 \text{ }^{\circ}\text{C}$

Deflasking time:  $\geq 90$  sec at  $37 \text{ }^{\circ}\text{C}$

Drewe Dentamid GmbH  
Max-Planck-Str. 31  
59423 Unna / Germany  
Tel.: 02303 / 8807-0  
Fax: 02303 / 8807-55  
dentamid@drewe.de  
www.drewe.com