

Technical Data Sheet

Product description: Bite registration material on vinyl polysiloxane basis.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	Chemical term
1	Polyvinylmethylsiloxane
2	Vinyl-Q Resin
3	Polymethyl hydrogensiloxane
4	Filling materials
5	Platinum catalyst
6	Pigments

Properties:

Colour: orange

Hardness: 48 ± 2 [Shore D] (DIN ISO 7619-1)

Consistency: Type 2: medium (DIN EN ISO 4823)

Recovery from deformation: > 99.9 % (DIN EN ISO 4823)

Linear dimensional change: < 0.1 % (DIN EN ISO 4823)

Processing time: 30 sec ± 5 sec at 23°C

Intraoral setting time: 80 sec ± 5 sec at 37°C
(incl. processing time)

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
dentamid@dreve.de
www.dreve.com